
136326-11

- Español
- Português
- Nederlands
- English

# QuikRead go<sup>®</sup> CRP

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





1. Price CP, Calvin J, Walker SA, Trull A, Newman DJ, Gorman EG. Rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. *Clin Chem Lab Med* 1999, 37, 109–113.
2. van Leeuwen MA, van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. *Baillière's Clinical Rheumatology* 1994, 531–552.
3. Olaison L, Hogevik H, Alestig K. Fever C-reactive protein and other acute-phase reactants during treatment of infective endocarditis. *Arch Intern Med* 1997, 157, 885–892.
4. Peltola H, Unkila-Kallio L, Kallio MJT. Simplified treatment of acute staphylococcal osteomyelitis of childhood. *Pediatrics*, 1997, 99, 846–850.
5. Philip AGS, Mills PC. Use of C-reactive protein in minimizing antibiotic exposure: Experience with infants initially admitted to a well-baby nursery. *Pediatrics*, 2000, 106, 1–5.
6. Pepys MB. The acute phase response and C-reactive protein. In: Warrell DA, Cox TM, Firth JD, Benz EJ, editors. *Oxford Textbook of Medicine*, 4th ed. Oxford University Press, 2003. Vol 2, p. 150–156.
7. Thomas L. C-reaktives Protein (CRP) In: *Labor und Diagnose*, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft, 1992, 780-790.

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

## Referenties • References

1. Price CP, Calvin J, Walker SA, Trull A, Newman DJ, Gorman EG. Rapid and sensitive automated light scattering immunoassay for serum C-reactive protein and the definition of a reference range in healthy blood donors. *Clin Chem Lab Med* 1999, 37, 109–113.
2. van Leeuwen MA, van Rijswijk MH. Acute phase proteins in the monitoring of inflammatory disorders. *Baillière's Clinical Rheumatology* 1994, 531–552.
3. Olaison L, Hogevik H, Alestig K. Fever C-reactive protein and other acute-phase reactants during treatment of infective endocarditis. *Arch Intern Med* 1997, 157, 885–892.
4. Peltola H, Unkila-Kallio L, Kallio MJT. Simplified treatment of acute staphylococcal osteomyelitis of childhood. *Pediatrics*, 1997, 99, 846–850.
5. Philip AGS, Mills PC. Use of C-reactive protein in minimizing antibiotic exposure: Experience with infants initially admitted to a well-baby nursery. *Pediatrics*, 2000, 106, 1–5.
6. Pepys MB. The acute phase response and C-reactive protein. In: Warrell DA, Cox TM, Firth JD, Benz EJ, editors. *Oxford Textbook of Medicine*, 4th ed. Oxford University Press, 2003. Vol 2, p. 150–156.
7. Thomas L. C-reactives Protein (CRP) In: *Labor und Diagnose*, 4th ed. Marburg: Die Medizinische Verlagsgesellschaft, 1992, 780-790.

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## 1 Intended use

The QuikRead go CRP test is intended for quantitative determination of CRP (C-reactive protein) in whole blood, serum and plasma using the QuikRead go instrument. For *in vitro* diagnostic use only.

## 2 Summary and explanation of the test

CRP is an acute phase protein present in low concentrations in healthy individuals<sup>1</sup>. Any pathological condition associated with invasive bacterial infection, inflammation or tissue destruction is accompanied by elevation of the CRP level in the patient's serum. The rise in CRP levels is rapid, and increased levels can be detected within 6 to 12 hours of the onset of the inflammatory process<sup>2</sup>.

Quantitative measurement of the CRP concentration has been reported to be a sensitive indicator of the efficacy of antimicrobial therapy and the course of bacterial infections, as well as an effective tool in controlling and monitoring postoperative infections<sup>2-6</sup>.

## 3 Principles of the procedure

QuikRead go CRP is an immunoturbidimetric test based on microparticles coated with anti-human CRP F(ab)<sub>2</sub> fragments. The CRP present in the sample reacts with the microparticles, and the resultant change in the turbidity of the solution is measured by the QuikRead go Instrument. QuikRead go CRP correlates well with results obtained using other immunoturbidimetric methods.

The sample is added to the buffer in prefilled cuvettes, and the blood cells in a whole blood sample are haemolysed. The calibration data of the test is in the bar code label of the cuvette, which QuikRead go reads automatically before the start of the test. The CRP value is automatically corrected according to the sample's hematocrit level. Hematocrit results are not displayed but used in calculations.

## 4 Reagents

### Content of the kit

Component name and origin	Symbol	QuikRead <sup>®</sup> go CRP Cat. no. 135171 50 tests
CRP Reagent Caps <b>ORIG</b> <b>SHP</b>	<b>REAG</b> <b>CPS</b>	2 x 25
Buffer in prefilled cuvettes	<b>BUF</b>	2 x 25 x 1 ml
Capillaries (20 µl)	<b>CAPIL</b> <b>HEP</b>	50
Plungers	<b>PLUN</b>	50
Instructions for use		

Reagents contain sodium azide, please see Section 5 "Warnings and precautions".

### Storage

Kit component	Storage at 2...8°C	Storage at 18...25°C
CRP Reagent caps (in opened and unopened aluminium tube)	Until the expiry date of the reagent caps	24 h per day – 1 month 7.5 h per day – 3 months
Prefilled cuvettes in unopened foil pouches	Until the expiry date of the prefilled cuvettes	Until the expiry date of the prefilled cuvettes
Unopened prefilled cuvettes without the foil pouch	6 months	3 months
Opened prefilled cuvette	2 hours	2 hours

Mark the date of opening the foil pouch on the cuvette rack.

### Reagent preparation and storage conditions

All reagents are ready to use. Keep the CRP Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.

### Reagent deterioration

The product should be used only if the buffer volume in the cuvette is correct. Check that the liquid surface is between the two lines marked on the cuvette. Do not use a cuvette with visible dirt in the buffer.

## 5 Warnings and precautions

### Health and safety information

- For *in vitro* diagnostic use only.
- Do not smoke, eat or drink in premises where samples or kit reagents are handled. Wear suitable personal protective clothing and disposable gloves when handling patient samples and kit reagents. Wash hands thoroughly when finished performing the test.
- Avoid contact with skin and eyes. After contact with skin, wash immediately with plenty of soap and water.
- All patient samples and controls should be handled as potentially infectious material.
- The buffer contains 0.004% reaction mass of: 5-chloro-2-methyl-2H-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-2H-isothiazol-3-one [EC no. 220-239-6] (3:1) (Skin Sens 1, Aquatic Chronic 3), and < 0.1% sodium azide. May cause an allergic skin reaction (H317). Harmful to aquatic life with long lasting effects (H412). Avoid breathing vapours (P261). Avoid release to the environment (P273). Wear protective gloves/protective clothing (P280). If skin irritation or rash occurs: Get medical advice/attention (P333+P313). Take off contaminated clothing and wash it before reuse (P362+P364). Dispose of contents according to national and local law (P501).
- Lyophilised reagents inside a reagent cap contain < 1% sodium azide (Aquatic Chronic 3). Harmful to aquatic life with long lasting effects (H412). Avoid release to the environment (P273). Dispose of contents according to national and local law (P501). Contact with acids liberates very toxic gas (EUH032).
- Reconstituted and liquid reagents contain < 0.1% sodium azide, which is not considered a harmful concentration. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Disposal: see Section 14.

**Warning**

### Analytical precautions

- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- The QuikRead go CRP kit is intended solely for use with the QuikRead go instrument.
- Do not mix components from different lot numbers or different tests. The components are disposable; never reuse components for performing a test.
- When opening a kit for the first time, ensure that the foil pouches protecting the cuvettes are intact. If a foil pouch is damaged, do not use the cuvettes inside it. Always ensure before using an individual cuvette that its foil covering is intact.
- Do not touch the clear flat surfaces on the lower part of the cuvette (optical part). Discard any cuvettes with fingerprints.
- QuikRead CRP Reagent Caps are colour-coded blue to distinguish them from other QuikRead analytes.

- Keep the QuikRead CRP Reagent Caps away from moisture. Close the aluminium tube immediately after taking out the required number of reagent caps.
- Do not splash liquid into the measuring well of the instrument.

## 6 Sample collection and preparation

### Sample material, collection and volume

Sample material	Sample volume	Sample collection
Fingertip whole blood	20 µl	Prick the clean and dry finger with a lancet and discard the first drop. Wipe the finger and collect 20 µl of blood from the second drop in a capillary.
Anticoagulated whole blood	20 µl	Whole blood collected in a tube containing heparin or EDTA can be used. Mix the whole blood by inverting the tube several times and take 20 µl with a capillary or pipette.
Plasma	20 µl or 12 µl *	Whole blood collected in a tube containing heparin or EDTA can be used. Avoid haemolysed samples. Mild haemolysis in samples will not affect the test results.
Serum	20 µl or 12 µl *	Avoid haemolysed samples. Mild haemolysis in samples will not affect the test results.

\*If you want to use sample volume 12 µl in order to have a wider measurement range, you have to change the settings of the QuikRead go instrument. Please see Instrument manual.

### Sample dilution

Sample material	Instructions
Whole blood	Do not dilute whole blood samples.
Plasma/serum	Plasma or serum samples may be diluted with a 0.9% NaCl solution before adding the samples to the cuvettes. Dilution ratio is 1+3 (1 part of volume sample + 3 parts of volume 0.9% NaCl). Take 20 µl of the diluted sample into a cuvette, analyse the sample and multiply the result by 4.

### Sample storage

Sample material	Short-term storage	Long-term storage
Fingertip whole blood	Not more than 15 minutes in the heparinised capillary	Do not store
Anticoagulated whole blood	2...8°C for 3 days	Separate plasma and store below -20°C
Plasma	2...8°C for 7 days	Below -20°C
Serum	2...8°C for 7 days	Below -20°C
Sample (whole blood, plasma, serum) in buffer	18...25°C for 2 hours	Do not store

Allow the samples to reach room temperature (18...25°C) before testing. Frozen samples must be thawed completely, mixed thoroughly and brought to room temperature prior to testing. If the thawed plasma samples contain clots the samples shall be centrifuged. Samples should not be repeatedly frozen and thawed.

## 7 Procedure

### Materials required but not supplied

Material	Cat. no.	
QuikRead go® Instrument	133893, 149915	
Capillaries (20 µl) 50 pcs	67962	Supplied with 135171
Plungers 50 pcs	67966	Supplied with 135171
QuikRead go® CRP Control	153764	Recommended for quality control
QuikRead go® CRP Control High	153763	Recommended for quality control
Fingertip lancets		

### Assay procedure

Discard the foil pouch of the cuvette rack and mark the opening date on the rack.

#### Sample collection (Fig. 1–5)

Glass capillaries (20 µl) and plungers are recommended for adding the sample into a prefilled cuvette. Capillaries (Cat. no. 67962) and plungers (Cat. no. 67966) can be purchased separately, whereas the QuikRead go CRP kit (Cat. no. 135171) contains them already.

- 1 A prefilled cuvette must reach room temperature (18...25°C) before use. This will take 15 minutes for an individual refrigerated (2...8°C) prefilled cuvette. Do not touch the clear flat surfaces in the lower part of the cuvette (the optical part). Remove the foil cover from the cuvette. Beware of splashing the liquid. Condensed buffer on the foil cover has no effect on results. The test must be performed within two (2) hours from opening the cuvette. The liquid surface should lie between the two lines marked on the cuvette.
- 2 Insert the plunger into the capillary at the end with the blue stripe.
- 3 Fill the capillary with the sample up to the white stopper. Ensure that there are no air bubbles in the capillary. Use a soft tissue to wipe away any residual sample from the outside of the capillary.
- 4 Place the sample capillary into the buffer solution in the cuvette and dispense by pressing down the plunger. Make sure that the capillary is completely empty.
- 5 Close the cuvette tightly with a CRP reagent cap. Do not press down the blue inner part of the cap. Once the sample is added to the buffer, it may be held for up to two (2) hours prior to running the assay. Keep the cuvette in an upright position, do not shake.

#### Analysing the sample (Fig. 6–8)

See the QuikRead go Instrument manual for detailed instructions. The display guides you in performing the test.

- 6 Choose **Measure** on the display of the QuikRead go instrument.
- 7 Insert the cuvette in the measurement well of the instrument. The barcode should be facing you as shown in the Figure 7. The instrument recognises the sample type. The display shows how the measurement is proceeding. First the instrument measures the sample blank. After that the instrument measures the CRP concentration within one minute.
- 8 When the measurement has been completed, the result will be shown in the display and the cuvette will rise up from the measurement well automatically.

#### Analysing the control (Fig. 7–8)

Please see the instructions for use provided with each control vial. Analyse control samples as patient samples but choose **Quality Control** on the display of the QuikRead go instrument (Figures 7–8). The result will be stored as the control measurement.

## 8 Quality control

Regular use of QuikRead go CRP Control (Cat. no. 153764) and/or QuikRead go CRP Control High (Cat. no. 153763) is recommended.

**QuikRead CRP Controls:** All QuikRead CRP controls are ready-to-use. The control value has been determined for the QuikRead go instrument with the same instructions for use as for plasma/serum samples. The sample volume is 20 µl. If you want to use sample volume 12 µl, you have to change the setting of the QuikRead go instrument. See Instrument manual. The control is stable for 15 minutes in the cuvette.

**Other commercially available CRP controls:** Handle and use according to the instructions for each control substance. Acceptable limits for the control must be defined with precision results gained with the QuikRead go instrument. The blanking process might not succeed if you use a control that contains artificial red blood cells.

## 9 Results' interpretation

Elevation of CRP value is unspecific and the results should be interpreted in the light of other clinical findings.

CRP test result	Interpretation of result <sup>7</sup>
< 10 mg/l	Excludes many acute inflammatory diseases, but does not specifically exclude inflammatory processes.
10–50 mg/l	Elevated concentrations found in acute disease that occur in the presence of a slight to moderate inflammatory process.
> 50 mg/l	Indicates high and extensive inflammatory activity.

## 10 Limitations of the procedure

Assay procedures other than those specified in these instructions may yield questionable results. Some substances may interfere with the test results; please see Section 12. "Performance Characteristics".

Test results should never be used alone, without a complete clinical evaluation, when making a diagnosis. Intra-individual variations in CRP are significant and should be taken into account – for example, by means of serial measurements – when interpreting the values.

## 11 Expected values

The assay reference range limit was determined by using 143 apparently healthy adults (59 males and 84 females) with age 19–65 according to CLSI EP28-A3C guideline. Based on the results the 95% reference limit was concluded to be  $\leq 5$  mg/l for each sample type. Results are based on a non parametric approach.

However, each laboratory is recommended to establish a range of normal values for their population in their region.

## 12 Performance characteristics

### Method comparison

Patient plasma samples were measured using two clinical laboratory methods, and the QuikRead go CRP method. Summary of correlation studies is presented in the table below.

Passing-Bablok analysis	
Clinical laboratory method 1	Clinical laboratory method 2
$y = 1.00x + 1.0$ $r = 0.994$ $n = 116$	$y = 0.93x + 1.4$ $r = 0.987$ $n = 113$

### Whole blood versus plasma comparability

In a comparison of 104 patient samples, whole blood (y) and plasma (x) were found comparable.  $y = 1.03x - 0.3$

### Precision

A precision study was performed according to Clinical and Laboratory Standards Institute (CLSI) guideline EP5-A2.

Within-run, between-day and total precision						
Sample material	Sample number	Number of days	Mean CRP (mg/l)	Within-run CV (%)	Between-day CV (%)	Total CV (%)
Whole blood	Sample 1	20	9	4.5	4.3	6.2
	Sample 2	20	52	1.4	2.4	3.2
	Sample 3	20	177	1.9	2.8	3.6
Plasma	Sample 1	20	10	3.7	0.1	4.2
	Sample 2	20	55	1.5	0.4	1.9
	Sample 3	20	148	2.0	1.0	2.4
Control	Sample 1	20	30	2.0	3.1	3.7

### Measuring range

For whole blood samples the measuring range is 5–200 mg/l at the normal hematocrit level of 40%. If the hematocrit is lower or higher than 40%, the measuring range will change according to the table below. Starting from QuikRead go instrument software version 3.1, the instrument displays a CRP result at hematocrit range of 15–75%. Previous instrument software versions display CRP results at hematocrit range of 20–60%. The CRP result is not displayed, if the level of hematocrit is outside the range of 15–75% or 20–60% (depending on the instrument software version).

Sample type	Hematocrit %	Sample volume	Measuring range mg/l CRP	
Whole blood sample	15–19	20 $\mu$ l	5–150	In cases when the result is above or below the measuring range, the result will be displayed, for example "> 200 mg/l CRP".
	20–28		5–160	
	29–35		5–180	
	36–41		5–200	
	42–46		5–220	
	47–50		5–240	
	51–53		5–260	
	54–56		5–280	
	57–60		5–300	
	61–62		5–330	
	63–64		5–340	
	65–66		5–360	
	67–68		5–390	
	69–70		5–410	
71–72	5–440			
73–74	5–470			
	75	5–510		
Plasma/ Serum sample	–	20 $\mu$ l	5–120	
		12 $\mu$ l	5–200	

### Interference

Interferent	No interference found up to concentration
Bilirubin	400 $\mu$ mol/l
Vitamin C	200 $\mu$ mol/l
Triglycerides	11.5 mmol/l
Cholesterol	9.0 mmol/l
Rheumatoid factors (RF)	525 IU/ml
Leucocytes	$145 \times 10^9$ cells
Anticoagulants (Li-heparin or EDTA)	No interference

Most heterophilic or anti-sheep antibodies in the samples do not interfere with the test, as the assay antibodies lack the Fc-part. In rare cases, interference by IgM myeloma protein has been observed.

### Antigen excess

CRP concentrations less than 1000 mg/l do not give falsely low results.

## 13 Traceability

The calibrators used to calibrate the CRP assay of the QuikRead go CRP test are traceable to the ERM®-DA474 reference material.

## 14 Disposal

- Dispose of contents according to national and local law.
- All patient samples, used caps, cuvettes, capillaries and plungers should be handled and disposed of as potentially infectious material.
- Materials of the components:  
 Paper: Instructions for use  
 Cardboard: Kit box  
 Plastic: Cuvettes, reagent caps, foil covering cuvette rack, cuvette rack, plungers, plunger and capillary tubes  
 Glass: Capillaries  
 Metal: Reagent cap tubes, cuvette lids, plunger and capillary tube caps  
 Several (not to be recycled): Lids of reagent cap tubes
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

## 15 Troubleshooting

### Error messages

Error messages of the QuikRead go instrument are listed in the table below. For more information about error messages, see the QuikRead go instrument manual.

Error messages	Corrective action
Cuvette position not correct. Remove cuvette.	Remove the cuvette and re-insert it into the right position.
Measurement prohibited. Please check the reagent cap.	Check that the cuvette has a reagent cap and the inner blue part of the cap is not pressed down.
Measurement prohibited. Kit lot expired.	Discard expired kit lot. Use a new one.
Measurement prohibited. Cuvette temperature too low.	Let the cuvette warm up to room temperature (18...25°C). Test the same cuvette again.
Measurement prohibited. Cuvette temperature too high.	Let the cuvette cool down to room temperature (18...25°C). Test the same cuvette again.
Test cancelled. Blank too high.	Test the same cuvette again. The blanking process has not been completed or the sample may contain interfering substances. In the latter case, the test cannot be completed.
Test cancelled. Unstable blank.	
Test cancelled. Error in reagent addition.	Perform a new test. There has been some problem during reagent addition. Make sure that the cap is properly closed.
Test cancelled. Instrument failure.	Perform a new test. If this message appears often, contact your local distributor or call customer service.
No result. Unable to do hematocrit correction.	Test the same cuvette again. If the error message appears again, perform a new test. If the test cannot be completed at all, the sample may contain interfering substances or be otherwise compromised.
No result. Hematocrit too low.	If the level of hematocrit in a whole blood sample is below the hematocrit measurement range of the instrument (see Section 12), no CRP result is shown.
No result. Hematocrit too high.	If the level of hematocrit in whole blood samples is above the hematocrit measurement range of the instrument (see Section 12), no CRP result is shown.
Error code xx. Please restart the instrument.	Restart the instrument. If the error message appears again, contact your local distributor or call customer service and give the error code.
Error code xx. Please call customer service.	Contact your local distributor or call customer service and give the error code.

### Unexpected low and high results

Possible reasons for unexpected low and high results are listed in the table below.

Problem	Possible cause	Corrective action
Unexpected low result.	Too small sample volume.	Perform a new test. Make sure that the capillary is completely filled. Avoid air bubbles.
Unexpected high result.	Too large sample volume.	Perform a new test. Make sure that the capillary is properly filled. Wipe away any excess sample from the outside of the capillary.
	Too low buffer volume.	Perform a new test. Make sure that the volume is correct by checking that the liquid surface is between the two lines marked on the cuvette.
	The cuvette is dirty.	Perform a new test. Do not touch the clear flat surfaces in the lower part of the cuvette.
Unexpected low/high result.	Component from different kit lots or tests are used.	Perform a new test. Make sure that all the reagents are from the same kit lot.
	Incorrect reagent storage.	Perform a new test. Make sure that the reagents are stored according to the instructions for use.
	Incorrect sample volume: plasma/serum/control sample	Check, that the used sample volume (20 or 12 µl) matches the set plasma/serum sample volume on the instrument.